

# LA-ICP-MS analysis of historic ink and paper

- ❑ **Soft laser ablation of historic documents**
- ❑ **Calibration standards for quantitative LA-ICP-MS analysis of ink and paper**

## Introduction

More than 200 years ago, the most commonly used writing ink in the Western World was an iron-gall ink, evident today in surviving manuscripts and drawings, including those by Leonardo da Vinci, Galileo, Rembrandt and many others.

Some of these historic documents are now suffering terrible damage, caused by breakdown of the paper, apparently accelerated by iron, impurity elements and other constituents in the ink. It has been found that Cu impurities in the iron sulfate, used in most ink preparations, is particularly damaging. To select the correct conservation strategy for such documents, it is essential to know the ink and paper compositions with respect to iron, copper and other transition elements.

## Aims of controlled low-power ablation

LA-ICP-MS yields sensitive elemental analysis. In the analysis of historic artifacts, discernible ablation marks on the analyzed area need to be avoided. Such "soft" ablation is achieved by the UP213 system, using optical attenuation of the beam output, as illustrated in Figure 2.

The aim of this work is to obtain quantitative LA-ICP-MS analysis of ink and paper by soft ablation sampling, requiring the preparation of synthetic paper analytical standards of known Group (II) and transition element content.

## Preparing the synthetic calibration standards

The synthetic paper standards, to be used for laser ablation calibrations, were prepared from purified cellulose paper (Whatman No. 1 filter paper, Maidstone, UK). Metals were introduced to the paper by immersion into a distilled water solution of 11 analytical grade metal chlorides (Mg, Ca, Ti, Mn, Fe, Co, Ni, Cu, Zn, Sr, and Cd).

The immersing solution of highest concentration contained 10 g/L of each metal salt and a sequence of diluted solutions was used for four other preparations. After a 2-minute immersion, each of the five paper samples was air-dried. The actual concentration of each metal in each paper was determined by solution ICP-MS analysis, following digestion.

## Analyzing the prepared synthetic standards

The final compositions of the paper standards were determined accurately after complete digestion of samples of each paper. Approximately 50 mg of each proposed paper standard was microwave-digested in  $\text{HNO}_3 / \text{H}_2\text{O}_2$ , yielding a clear, particle free solution.<sup>(1)</sup> Each of the five papers was sampled and analyzed in triplicate and the same digestion procedure was also used on five replicates of pure Whatman paper, used for blank determination.



Fig 1. Many historic documents contain iron-gall ink




Fig 2. "Soft" laser ablation line of the UP 213 laser.

## Instrumentation

New Wave UP213 laser system

Ablation gas: He at 0.95 L/min  
Plasma make-up gas: Ar at 0.75 l/min  
Crater diameter: 100 µm  
Scan speed: 25 µm/s  
Repetition rate: 10 Hz  
Energy density: 0.1 J/cm<sup>2</sup>

## LA-ICP-MS calibration

Having determined the composition of all five standards, they were then used to calibrate LA-ICP-MS response. Each standard (5 mm x 7 mm) was mounted onto a microscopic glass slide with self-adhesive tape. Six consecutive replicate acquisitions of each calibration sample and the blank were measured and  $^{13}\text{C}$  was used as an internal standard.

LA-ICP-MS calibrations show that there is excellent correlation between LA-ICP-MS signal and the composition of the standard paper. All 11 elements yielded linear calibration plots of this general pattern.




Fig 3. Calibration graphs for Ca, Mg, Zn and Fe.

## Analysis of an unknown sample

The calibration standards were then used in LA-ICP-MS analysis of an 18<sup>th</sup> Century document and of a laboratory model document. The model document's ink was made using a known 18<sup>th</sup> Century formulation.

The laboratory model paper was the same cellulose paper used to prepare the standard papers, with no sizing used. The analysis summary, in Table 1, shows that the historic ink contains significant minor element impurities that do not occur in the newly formulated ink. The historic ink was only required to have good colour and adhesion to the paper. Ingredients available at that time were clearly not high in purity.

Element content [mg/g]	Original document from 18th century		Laboratory model document	
	paper	ink line	paper	ink line
Mg	0.61 ± 0.09	<b>4.3 ± 0.7</b>	0.13 ± 0.05	0.13 ± 0.02
Ca	0.7 ± 0.2	<b>5.0 ± 2</b>	0.46 ± 0.03	0.44 ± 0.04
Ti	< 0.05	< 0.05	< 0.05	< 0.05
Mn	0.21 ± 0.06	<b>0.4 ± 0.1</b>	0.18 ± 0.02	0.32 ± 0.05
Fe	0.30 ± 0.09	<b>39 ± 10</b>	0.23 ± 0.06	<b>1.4 ± 0.2</b>
Co	0.19 ± 0.06	0.20 ± 0.08	0.19 ± 0.08	0.19 ± 0.03
Ni	0.3 ± 0.1	0.3 ± 0.1	0.3 ± 0.1	0.3 ± 0.1
Cu	0.3 ± 0.1	0.35 ± 0.09	0.26 ± 0.03	0.36 ± 0.05
Zn	0.08 ± 0.04	<b>7 ± 2</b>	0.05 ± 0.02	0.07 ± 0.03
Sr	0.24 ± 0.06	0.28 ± 0.09	0.24 ± 0.07	0.24 ± 0.03
Cd	< 0.05	< 0.05	< 0.05	< 0.05

Table 1: Comparison of original document and laboratory model. Significant levels in inks are shown in bold.




Fig 4. Laser traverse along ink line

Figure 4 shows Fe, Zn and Ca signals during a traverse across an ink line. Even though the UP 213 generates a barely discernible ablation track, the minor element content of the ink can be clearly differentiated from that of the paper.

## Conclusions

Excellent LA-ICP-MS calibrations have been achieved for 11 elements using synthetic paper calibration standards.

An 18<sup>th</sup> Century iron-based ink has been shown to contain significant levels of Zn, Mn, Ca and Mg. A modern ink, made from higher purity compounds, has fewer such concomitant elements.

Using very low laser energy density, ablations tracks were only visible using a microscope.

## Acknowledgment

Vid Simon Šelih  
University of Ljubljana, Faculty of Chemistry and Chemical Technology, Ljubljana, Slovenia.

## References

(1) Milestone application note DG-CH-07.



**USA**  
New Wave Research, Inc  
48660 Kato Road  
Fremont CA 94538-7339  
Tel: 510-249-1550  
Tel: 800-566-1743  
Fax: 510-249-1551  
Email: [Lasers@new-wave.com](mailto:Lasers@new-wave.com)

**Japan**  
New Wave Research, KK  
5F Chojiya Building, 1-36-4,  
Shinjuku-ku, Shinjuku  
Tokyo, 160-0022 Japan  
Tel: +81-3-3351-0131  
Fax: +81-3-3351-0121  
Email: [NewWaveKK@new-wave.com](mailto:NewWaveKK@new-wave.com)

**Taiwan**  
New Wave Research G. C. Co., Ltd.  
2Fl., No. 58-3, Xing Shan Road,  
Neihu Dist., Taipei  
Taiwan 114  
Tel: 886-2-8792-7585  
Fax: 886-2-8792-7584  
Email: [NewWaveCG@new-wave.com](mailto:NewWaveCG@new-wave.com)

**Europe**  
New Wave Research Co. Ltd.  
8 Avro Court  
Ermine Business Park  
Huntingdon, Cambridge  
PE29 6XS, UK  
Tel: 44-(0)1480 456 566  
Fax: 44-(0)1480 456 545  
Email: [NewWaveEU@new-wave.com](mailto:NewWaveEU@new-wave.com)

**China**  
New Wave Research (China)  
Co., Ltd  
Rm. 1102, 11/F, Jiaying  
Mansion,  
No.877, Dongfang Rd.  
Shanghai 200122, China  
Tel: 86-21-5058-7785  
Fax: 86-21- 5058-7786  
Email: [NewWaveCN@new-wave.com](mailto:NewWaveCN@new-wave.com)